

	<b>BIO-PM 4.14 Thermometer Calibration</b>	
	<i>Document #: 1921</i>	<i>Page 1 of 2</i>
	<i>Revision #: 2</i>	<i>Issued Date: 12/12/2018</i>
	<i>Document Manager: Kristin Schelling</i>	<i>Approved By: Jeffrey Nye</i>

## 4.14 Thermometer Calibration

### 4.14.1 Scope

To accurately calibrate thermometers used in the DNA laboratory.

### 4.14.2 Background Information

NIST Manual

### 4.14.3 Operation

#### 4.14.3.1

The report on calibration of the NIST thermometer shows that accurate thermometer readings are at 3 basic temperatures:

Reading	Correction
0.08C	-0.08C
37.00	0.10
56.00	0.15

#### 4.14.3.2

The comparison of these temperatures of the NIST thermometer to a second thermometer will allow for calibration of that second thermometer.

#### 4.14.3.3

Calibration of the thermometers or digital readouts for the instruments will be recorded on designated forms.

### 4.14.4 Maintenance

#### 4.14.4.1

Store thermometer in the case provided. Maintain the NIST certification information provided when thermometers are purchased.

#### 4.14.4.2

Use caution in prevention of the thermometer mercury from separating.

	<b>BIO-PM 4.14 Thermometer Calibration</b>	
	<i>Document #: 1921</i>	<i>Page 2 of 2</i>
	<i>Revision #: 2</i>	<i>Issued Date: 12/12/2018</i>
	<i>Document Manager: Kristin Schelling</i>	<i>Approved By: Jeffrey Nye</i>

#### **4.14.4.3**

In case of separation of the mercury, refer to the operator manuals for recombination techniques.