

	BIO-PM 4.34 Dri-Bath and Heat Block	
	<i>Document #: 1942</i>	<i>Page 1 of 1</i>
	<i>Revision #: 1</i>	<i>Issued Date: 11/14/2017</i>
	<i>Document Manager: Kristin Schelling</i>	<i>Approved By: Jeffrey Nye</i>

4.34 Dri-Bath And Heat Block

4.34.1 Scope

To provide a controlled dry heat environment for specific sized test tubes and/or microplates.

4.34.2 Background Information

See Thermolyne Dri-Bath operator's manual.

4.34.3 Operation

4.34.3.1

Insert a 45 mm immersion depth thermometer into the block.

4.34.3.2

Turn control switch ON.

4.34.3.3

When temperature level of thermometer is on the desired setting, slight adjustments right or left after stabilization of the temperature will produce the correct temperature.

4.34.4 Maintenance

4.34.4.1

Clean both the Dri-bath holder and the aluminum blocks with a mild soap and water solution.

4.34.4.2

Should control of the heating mechanism fail-replace the unit.