

	ARSON-PM 5.1.6 Fire Debris Analysis - Sample Extract Storage	
	Document #: 3115	Page 1 of 1
	Revision #: 1	Issued Date: 12/28/2017
	Document Manager: Cheryl Lozen	Approved By: Jeffrey Nye

A sampling of the charcoal strip extracts from fire debris samples and questioned liquids analyzed for the presence of ignitable liquid residues will be preserved for potential reanalysis. An extract sampling will be returned to the agency with the evidence. See 5.1.2.3 Static (Passive) Headspace Concentration (Adsorption-Elution) - **Extract Handling** for instructions.

Extract Storage Short-Term---Extracts, liquid samples, or solvents may be stored in septum seal or screw cap glass vials with PTFE-lined seals to prevent evaporation.